



# ***SHAKER / HOMOGENEIZER IN 'V'***

## ***TE-201/05***

Used to homogenize samples in powder or granules.

## Technical Characteristics

### TE-201/05

- Agitation range: 360°;
- Homogenizer: Built in 304 stainless steel with polished finish, internal and external;
- Gearmotor: Induction with fixed rotation of  $28 \pm 2$  RPM;
- Covers: In polished stainless steel 304, removable, with silicone seal;
- Capacity: 15 Liters;
- Cabinet: In carbon steel with anti-corrosive treatment and electrostatic painting;
- Overall dimensions: W = 940 x D = 420 x H = 850 mm;
- Weight: 44 kg;
- Power: 245 Watts;
- Voltage: 220 Volts;
- Accompanies: - 02 extra fuses - Instruction Manual with Warranty Term;
- Note: The reported volume is based on a product with density 1;
- Comments: See other models with different volumes;

## Benefits and Advantages

- Easy-to-remove lid at the three ends for cleaning, sample entry and material exit
- 304 stainless steel bowl with solder and fine polishing to ensure a good asepsis, avoiding contamination and for a longer useful life
- Brushless induction motor, which can be used for long periods
- Strict Quality Control, in which checks and tests guarantee the perfect functioning of the equipment, providing safety and customer satisfaction
- Customer service, to answer questions and provide explanations about the equipment and methodologies
- Possibility of adaptations according to the needs of the customer, makes the equipment already in line a special equipment.