



PIPETTOR SEMI AUTOMATIC FOR SOIL TECSOLO-005

Used for analysis in soil preparation for fertility purposes.

Technical Characteristics

TECSOLO-005

- Capillary pipettes: 11 borosilicate pipettes;
- Volume: 5 ml;
- Structure: fully stainless steel;
- Dimensions: L = 918 X D=320 x H = 810 mm;
- Weight: 16kg;
- Accompany: - 45 Clamps Warranty Term NOTE: it is indispensable to use vacuum pump;

Benefits and Advantages