



SEMI-AUTOMATIC PIPETTE

TECSOLO-050

Used for analysis in soil preparation, for fertility purposes.

Technical Characteristics

TECSOLO-050

- Capillary pipettes: 11 borosilicate pipettes;
- Volume: 50 ml;
- Structure: I fully stainless steel;
- Dimensions: L = 918 X D=320 x H = 910 mm;
- Weight: 16kg;
- ACCOMPANY: - 45 Clamps - Warranty term;
- Note: it is indispensable to use vacuum pump;

Benefits and Advantages