



# ***PIPETTOR SEMI AUTOMATIC***

## ***TECSOLO-100***

Used for analysis in soil preparation for fertility purposes.

## Technical Characteristics

### TECSOLO-100

- Capillary pipettes: 11 borosilicate pipettes;
- Volume: 100ml;
- Structure: fully stainless steel;
- Dimensions: L = 918 X D=320 x H=1015 mm;
- Weight: 16kg;
- ACCOMPANY: - 45 Clamps - Instruction Manual with Warranty Term;
- Note: it is indispensable to use vacuum pump;

Benefits and Advantages