



# ***TITRATION BRIDGE***

## ***TECSOLO-200***

Used to aid solution titrations in soil analysis.

## Technical Characteristics

### TECSOLO-200

- Capacity: 11 125ml erlenmeyers or pH cups;
- Support structure: fully stainless steel;
- Base and sliding structure: All in stainless steel;
- Dimensions: L = 960 x d=400 x H=965 mm;
- Weight: 10kg with shaker;
- Power: 50 Watts;
- Voltage: 220 Volts;
- ACCOMPANY: - 01 borosilicate glass Burette with capacity for 10 mL - 01 Magnetic shaker without heating - Warranty term;

### Benefits and Advantages